

Implant Model

Detailed planning models for implant cases.

Print solid, accurate models with soft-tissue and analog options - ideal for prosthetic planning and patient communication.

Workflow Timeline

Total Time: ~50-70 min

1	Scan Capture arch and soft tissue with ODYX S1.	~1 min
2	Design Prepare model bases and analog sockets in CAD.	~6 min
3	Print Print models on P1-26 with Model Resin.	~35 min
4	Cure Wash and cure for dimensional stability.	~8 min
5	Assemble Seat analogs and verify soft-tissue fit.	~5 min

Printing Parameters (P1-26)

Layer Thickness	50 um
Exposure Time	1.8 - 2.5 s
Bottom Exposure	20 - 28 s
Bottom Layers	4 - 5
Build Plate Temp.	28 - 32 degC
Resin	Model Resin

Clinical Tips

- Orient models to protect critical soft-tissue areas.
- Verify analog torque seating before wax-up or CAD.
- Label arches clearly for multi-unit cases.
- Store cured models away from direct light.