

Retainers

Retention appliances, printed on demand.

Produce accurate models for clear retainers and retention appliances - fast replacements without new impressions.

Workflow Timeline

Total Time: ~50-60 min

1	Scan Capture final occlusion with ODYX S1.	~45 sec
2	Design Prepare retainer model bases in CAD.	~3 min
3	Print Print models on P1-26 with Model Resin.	~30 min
4	Cure Wash and cure before thermoforming.	~8 min
5	Thermoform Form and trim the retainer.	~12 min

Printing Parameters (P1-26)

Layer Thickness	50 um
Exposure Time	2.5 - 3.0 s
Bottom Exposure	25 - 30 s
Bottom Layers	5
Build Plate Temp.	30 - 35 degC
Resin	Model Resin

Clinical Tips

- Archive final scans for lifetime retainer remakes.
- Check occlusion contacts after thermoforming.
- Polish edges for patient comfort.
- Provide care instructions with every delivery.