

Surgical Guide

Accurate implant placement, digitally planned.

Print rigid, precise surgical guides from CBCT-driven plans - transfer the digital plan to the chair with confidence.

Workflow Timeline

Total Time: ~45-60 min

1	Scan and CBCT Merge IOS and CBCT for implant planning.	~2 min
2	Design Plan implants and design the guide sleeves.	~10 min
3	Print Print guide on P1-26 with Surgical Guide Resin.	~25 min
4	Cure Wash and fully cure for clinical rigidity.	~8 min
5	Verify Check fit on model, then sterilize per protocol.	~5 min

Printing Parameters (P1-26)

Layer Thickness	50 - 100 um
Exposure Time	2.0 - 2.8 s
Bottom Exposure	25 - 35 s
Bottom Layers	5 - 6
Build Plate Temp.	30 - 35 degC
Resin	Surgical Guide Resin

Clinical Tips

- Validate guide fit on the printed model before surgery.
- Confirm sleeve diameter against the surgical kit.
- Follow wash/cure protocols for biocompatible resins.
- Inspect for supports residue near sleeve openings.