

Try-ins

Perfect try-ins for better fit and function.

Print inexpensive try-in appliances to validate esthetics, occlusion, and phonetics before committing to the final prosthesis.

Workflow Timeline

Total Time: ~60-75 min

1	Scan Record arches and occlusal relationship.	~2 min
2	Design Design try-in from the proposed setup.	~12 min
3	Print Print try-in on P1-26 with Temporary Resin.	~30 min
4	Cure Cure enough for clinical try-in strength.	~8 min
5	Evaluate Assess fit, esthetics, and phonetics.	~10 min

Printing Parameters (P1-26)

Layer Thickness	50 um
Exposure Time	2.5 - 3.0 s
Bottom Exposure	25 - 30 s
Bottom Layers	5
Build Plate Temp.	30 - 35 degC
Resin	Temporary Resin

Clinical Tips

- Photograph try-in for lab or in-house design notes.
- Check midline, plane, and lip support carefully.
- Mark adjustments on the try-in before redesign.
- Do not use try-in resin as a long-term prosthesis.