

Aligners

Clear aligner workflows, end to end.

Print precise thermoforming models for clear aligner series - consistent staging, efficient batch production.

Workflow Timeline

Total Time: per series batch

1	Scan Full-arch scan with ODYX S1 for staging.	~1 min
2	Setup Plan movements in your aligner software.	~15 min
3	Print Batch-print staging models on P1-26.	~40 min
4	Cure Wash and cure models before thermoforming.	~10 min
5	Thermoform Form, trim, and polish aligners.	~20 min

Printing Parameters (P1-26)

Layer Thickness	50 - 100 um
Exposure Time	1.8 - 2.4 s
Bottom Exposure	20 - 28 s
Bottom Layers	4 - 5
Build Plate Temp.	28 - 32 degC
Resin	Model Resin

Clinical Tips

- Nest models efficiently to maximize build plate use.
- Fully dry models before thermoforming sheets.
- Inspect attachments on each staging model.
- Label stages clearly to avoid sequence errors.