

Dentures

Complete and partial dentures, digitally.

Print bases and try-in components for efficient full and partial denture workflows - predictable fit, fewer appointments.

Workflow Timeline

Total Time: depends on case type

1	Scan Capture arches and bite with ODYX S1.	~2 min
2	Design Design base and tooth setup in CAD.	~20 min
3	Print Print bases/teeth on P1-26.	~45 min
4	Cure Wash and cure for biocompatibility.	~12 min
5	Finish Polish, characterize, and deliver.	~15 min

Printing Parameters (P1-26)

Layer Thickness	50 - 100 um
Exposure Time	2.2 - 3.2 s
Bottom Exposure	28 - 40 s
Bottom Layers	5 - 6
Build Plate Temp.	30 - 35 degC
Resin	Denture Base Resin

Clinical Tips

- Verify vertical dimension at try-in before final cure finish.
- Follow manufacturer wash/cure for mucosal contact resins.
- Check border extension for comfort and retention.
- Document shade and tooth mold for remakes.