

Inlays and Onlays

Precise partial restorations, chairside.

Restore cusps and proximal anatomy with accurately printed inlays and onlays - conservative prep, excellent margins.

Workflow Timeline

Total Time: ~30-45 min

1	Scan Quadrant scan with clear prep margins.	~30 sec
2	Design Design inlay/onlay contacts and occlusion.	~4 min
3	Print Print on P1-26 at recommended layer height.	~16 min
4	Cure Fully cure before try-in and cementation.	~5 min
5	Finish and Cement Adjust contacts lightly, then cement.	~6 min

Printing Parameters (P1-26)

Layer Thickness	50 um
Exposure Time	2.5 - 3.0 s
Bottom Exposure	25 - 30 s
Bottom Layers	5
Build Plate Temp.	30 - 35 degC
Resin	Ceramic Crown Resin

Clinical Tips

- Capture sharp prep margins with powder-free scanning.
- Check proximal contacts on a solid model before cementation.
- Avoid over-thinning occlusal tables beyond material limits.
- Use dual-cure resin cement for deep preparations.