

Same-Day Crown

From scan to cementation in one visit.

Deliver strong, esthetic crowns in a single appointment with the ODYX digital workflow.

Workflow Timeline

Total Time: ~30-45 min

1	Scan Full arch or quadrant scan with ODYX S1.	~30 sec
2	Design Design the crown with your preferred CAD software.	~3 min
3	Print Print the crown with ODYX P1-26 using Ceramic Crown Resin.	~18 min
4	Cure Cure with ODYX Cure for optimal strength and esthetics.	~5 min
5	Finish and Cement Simple finishing and polishing. Ready to cement.	~5 min

Printing Parameters (P1-26)

Layer Thickness	50 um
Exposure Time	2.5 - 3.0 s
Bottom Exposure	25 - 30 s
Bottom Layers	5
Build Plate Temp.	30 - 35 degC
Resin	Ceramic Crown Resin

Clinical Tips

- Ensure proper isolation and shade selection.
- Minimal adjustments needed after try-in.
- Use fine diamonds and polishing kit for best results.
- Follow recommended curing time for maximum properties.