

Veneers

Esthetic veneers in a single visit.

Print thin, highly esthetic veneers chairside with the ODYX digital workflow - precise fit, natural translucency.

Workflow Timeline

Total Time: ~35-50 min

1	Scan Capture prep and antagonist with ODYX S1.	~45 sec
2	Design Design veneer morphology and emergence in CAD.	~5 min
3	Print Print veneers on P1-26 with Ceramic Crown Resin.	~15 min
4	Cure Post-cure for strength and optical properties.	~5 min
5	Finish and Bond Polish and bond with conventional protocols.	~8 min

Printing Parameters (P1-26)

Layer Thickness	50 um
Exposure Time	2.5 - 3.0 s
Bottom Exposure	25 - 30 s
Bottom Layers	5
Build Plate Temp.	30 - 35 degC
Resin	Ceramic Crown Resin

Clinical Tips

- Verify shade under multiple light sources before bonding.
- Keep veneer thickness within material indications.
- Use a calibrated try-in paste for shade confirmation.
- Polish margins carefully to preserve enamel seal.